

Product Datasheet

EAAC1 antibody [Out of stock] (orb149920)

Description

Mouse monoclonal antibody against

Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	ATTO-488
Tested Applications	ICC, IP, WB
Immunogen	Fusion protein amino acids 1-524 (full length) of human EAAC1
Target	EAAT3
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	5 µg of SMC-406 was sufficient for immunoprecipitation of EAAC1 from 50 µg of rat brain lysate.
Isotype	IgG1
Clonality	Monoclonal
MW	60kDa
Uniprot ID	P43005
NCBI	NP_004161.4
Entrez	6505
Dilution Range	WB (1:1000), ICC/IF (1:100);optimal dilutions for assays should be determined by the user.
Expiration Date	12 months from date of receipt.