

Product Datasheet

**Human CXCL2 Protein, Fc Tag (MALS verified)
(orb1499260)**

Description

Human CXCL2, Fc Tag is expressed

contains AA Ala 35 - Asn 107 (Acc...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

CXCL2

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Tag

C-Fc

Note

For research use only

Application notes

This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 34.3 kDa. The protein migrates as 38-42 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Purity

95%

Source

Human CXCL2, Fc Tag (orb1499260) is expressed from human 293 cells (HEK293). It contains AA Ala 35 - Asn 107 (Accession # P19875).

MW

34.3 kDa

Uniprot ID
P19875
Expiration Date

6 months from date of receipt.