

Product Datasheet

Alpha Tubulin (Acetyl-Lys40) antibody (orb1499276)

Description	Alpha Tubulin (Acetyl-Lys40) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Alpha Tubulin (Acetyl-Lys40) (SD(Ac-K)TI)
Target	Alpha Tubulin (Acetyl-Lys40)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	50 kDa
Uniprot ID	P68366
Expiration Date	12 months from date of receipt.