

Product Datasheet

AKT1/AKT2/AKT3 antibody (orb1499280)

Description	AKT1/AKT2/AKT3 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human AKT1 / AKT2 / AKT3
Target	AKT1 / AKT2 / AKT3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	56 kDa
Uniprot ID	P31751
Expiration Date	12 months from date of receipt.