

Product Datasheet

MYH7B antibody (orb1499345)

Description	MYH7B antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein within C-terminal Human MYH7B.
Target	MYH7B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	226 kDa
Uniprot ID	A7E2Y1
Expiration Date	12 months from date of receipt.