

Product Datasheet

TAK1/MAP3K7 antibody (orb1499366)

Description	TAK1/MAP3K7 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	KLH conjugated synthetic peptide derived from human TAK1/MAP3K7
Target	TAK1/MAP3K7
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	67 kDa
Uniprot ID	O43318
Expiration Date	12 months from date of receipt.