

Product Datasheet

GATA1 antibody (orb1499369)

Description	GATA1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant protein within human GATA1 (1-250/413 aa)
Target	GATA1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	45 kDa
Uniprot ID	P15976
Expiration Date	12 months from date of receipt.