

Product Datasheet

beta Tubulin antibody (orb1499401)

Description	beta Tubulin antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human beta Tubulin
Target	beta Tubulin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Recombinant
MW	50 kDa
Uniprot ID	P07437
Dilution Range	WB=1:500-1000, IHC-P=1:100-500, ICC=1:20-100, IF=1:20-100
Expiration Date	12 months from date of receipt.