

Product Datasheet

DNA Polymerase beta antibody (orb1499417)

Description	DNA Polymerase beta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant human DNA Polymerase beta protein (1-335/335 aa)
Target	DNA Polymerase beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P06746
Dilution Range	WB=1:200-1000
Expiration Date	12 months from date of receipt.