

Product Datasheet

Clcn3 antibody (orb149967)

Description

Mouse monoclonal to Clcn3 (PerCP). Clcn3

protons, and...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

PerCP

Tested

ICC, IHC, WB

Applications

Immunogen

Synthetic peptide amino acids 98-115
(cytoplasmic N-terminus) of rat Clcn3

Target

CICN3

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-408 was sufficient for detection of Clcn3 in 20 µg of rat brain membrane lysate and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

90kDa

Uniprot ID

[P51792](#)

NCBI

[NP_445815.2](#)

Entrez

[84360](#)

Dilution Range

WB (1:1000), IHC (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.