

Product Datasheet

SUR1 antibody (orb149980)

Description

Mouse monoclonal to SUR1 (Biotin). Sulfonylurea

	are the ...
Species/Host	Mouse
Reactivity	Hamster, Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ICC, IHC, WB
Immunogen	Fusion protein amino acids 1548-1582 (cytoplasmic C-terminus) of rat SUR1
Target	SUR1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-409 was sufficient for detection of SUR1 in 20 µg of mouse brain membrane lysate and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	160kDa
Uniprot ID	Q09429
NCBI	NP_037171.2
Entrez	25559
Dilution Range	WB (1:1000), IHC (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.