



Product Datasheet

eIF4G (phospho-Ser1231) Antibody (orb15002)

Description	eIF4G (phospho-Ser1231) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q04637
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of eIF4G(phospho-S...