

Product Datasheet

GIT1 antibody (orb150047)

Description

Mouse monoclonal to GIT1 (APC). G-protein

proteins 1 and ...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

APC

Tested

ICC, IP, WB

Applications
Immunogen

Fusion protein amino acids 375-770 (C-terminus)
of rat GIT1

Target

GIT1

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-413 was sufficient for detection
of GIT1 in 10 µg of rat brain lysate by
colorimetric immunoblot analysis using Goat
anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

90kDa

Uniprot ID

[Q9Z272](#)

NCBI

[NP_114002.1](#)

Entrez

[83709](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.