

Product Datasheet

Dicer antibody [Out of stock] (orb150089)

Description

Mouse monoclonal to Dicer (ATTO-390).

that specifically cl...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-390

Tested

ICC, WB

Applications
Immunogen

Fusion protein amino acids 1638-1899 of mouse Endoribonuclease Dicer

Target

Dicer

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-416 was sufficient for detection of Dicer in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

215kDa

Uniprot ID

[Q8R418](#)

NCBI

[NP_683750.2](#)

Entrez

[192119](#)

Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.