



Product Datasheet

c-kit (phospho-Tyr936) Antibody (orb15011)

Description	c-kit (phospho-Tyr936) Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P10721
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Western blot analysis of c-Kit(phospho-T...