

Product Datasheet

**LRP4 (Extracellular) antibody [Out of stock]
(orb150126)**

Description

Mouse monoclonal to LRP4 (ATTO-

receptor, which binds and inte...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-594

Tested

ICC, WB

Applications
Immunogen

Fusion protein amino acids 26- 350 (extracellular N-terminus) of mouse LTP4

Target

LRP4

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-418 was sufficient for detection of LRP4 in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2A

Clonality

Monoclonal

MW

215kDa

Uniprot ID

[Q8VI56](#)

NCBI

[NP_766256.3](#)

Entrez

[83469](#)

Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.