

Product Datasheet

TARP Gamma2/4/8 (Stargazin) antibody [Out of stock] (orb150192)

Description

Mouse monoclonal to TARP

belongs to a group called voltage gated calciu...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-488

Tested

ICC, WB

Applications
Immunogen

Fusion protein amino acids 203-323 (Cytoplasmic C-terminus) of rat TARPGamma2

Target

Stargazin

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-422 was sufficient for detection of TARP Gamma2/4/8 in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2A

Clonality

Monoclonal

MW

35-55kDa

Uniprot ID
[Q71RJ2](#)
NCBI
[NP_445803.1](#)
Entrez
[84347](#)
Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.