

Product Datasheet

TARP Gamma2/4/8 (Stargazin) antibody [Out of stock] (orb150199)

Description

Mouse monoclonal to TARP

	Stargazin belongs to a group called voltage ...
Species/Host	Mouse
Conjugation	AP
Tested Applications	ICC, WB
Immunogen	Fusion protein amino acids 203-323 (Cytoplasmic C-terminus) of rat TARP γ 2
Target	Stargazin
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 μ g/ml of SMC-422 was sufficient for detection of TARP γ 2/4/8 in 20 μ g of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2A
Clonality	Monoclonal
MW	35-55kDa
Uniprot ID	Q71RJ2
NCBI	NP_445803.1
Entrez	84347
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.