



Product Datasheet

JNK1/JNK2/JNK3 (phospho-Thr183/Tyr185) Antibody (orb15028)

Description	JNK1/JNK2/JNK3(phospho-Thr183/Tyr185) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P45983
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescent
analysis of H...



Western blot
analysis of
JNK1/JNK2/JNK3(...)