


## Product Datasheet

# EGFR (phospho-Tyr1197) Goat Polyclonal Antibody (orb15029)

|                            |                                                                                           |
|----------------------------|-------------------------------------------------------------------------------------------|
| <b>Description</b>         | EGFR (phospho-Tyr1197) Goat Polyclonal Antibody                                           |
| <b>Species/Host</b>        | Goat                                                                                      |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                              |
| <b>Tested Applications</b> | WB                                                                                        |
| <b>Storage</b>             | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| <b>Note</b>                | For research use only                                                                     |
| <b>Clonality</b>           | Polyclonal                                                                                |
| <b>Uniprot ID</b>          | <b>P00533</b>                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                           |



Western blot analysis of EGFR(phospho-Ty...