

Product Datasheet

ASIC1 antibody (orb150291)

Description

Mouse monoclonal to ASIC1 (RPE). Acid sensing

	BNAC2 and AS...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	RPE
Tested Applications	ICC, IHC, WB
Immunogen	Fusion protein amino acids 460-526 (Cytoplasmic C-terminus) of mouse ASIC1
Target	ASIC1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-427 was sufficient for detection of ASIC1 in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	60kDa
Uniprot ID	Q6NXX8
NCBI	NP_033727.1
Entrez	11419
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.