

Product Datasheet

Brevican antibody (orb150303)

Description

Mouse monoclonal to Brevican (Biotin).

sulfate proteoglyca...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Biotin

Tested

ICC, WB

Applications

Immunogen

Fusion protein amino acids 219-655 of rat Brevican

Target

Brevican

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-428 was sufficient for detection of Brevican in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2B

Clonality

Monoclonal

MW

140kDa

Uniprot ID

[P55068](#)

NCBI

[NP_001028837.1](#)

Entrez

[25393](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.