

Product Datasheet

SUR2A antibody (orb150345)

Description

Mouse monoclonal to SUR2A (ATTO-488). Sulfonylurea

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-488
Tested Applications	ICC, IHC, WB
Immunogen	Fusion protein amino acids 1505-1546 (SSIVDAGLVLVFSEGLVECDTGPNLLQHKNGLFSTLVMTNK, cytoplasmic C-terminus) of mouse SUR2A
Target	SUR2A
Preservatives	PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-431 was sufficient for detection of SUR2A in 20 µg of mouse brain membrane lysate and assayed by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2A
Clonality	Monoclonal
MW	120kDa
Uniprot ID	P70170
NCBI	NP_001038185.1
Entrez	20928
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.