

Product Datasheet

Versican antibody (orb150495)

Description

Mouse monoclonal to Versican (RPE). Versican

extrace...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

RPE

Tested

ICC, WB

Applications

Immunogen

Fusion protein amino acids 362-585
(glycosaminoglycan alpha domain) of mouse
Versican core protein

Target

Versican

Preservatives

PBS pH 7.4, 50% glycerol, 0.1% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-439 was sufficient for detection
of Versican in 20 µg of mouse brain membrane
lysate and assayed by colorimetric immunoblot
analysis using goat anti-mouse IgG:HRP as the
secondary antibody.

Isotype

IgG1

Clonality

Monoclonal

MW

350kDa

Uniprot ID

[Q62059](#)

NCBI

[AAH96495](#)

Entrez

[13003](#)

Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.