

Product Datasheet

Ataxin 1 antibody [Out of stock] (orb150770)

Description

Mouse monoclonal to Ataxin 1 (ATTO-594).

	family and contai...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-594
Tested Applications	ICC, IHC, IP, WB
Immunogen	Synthetic peptide amino acids 164-197 (ATTPSQRSQLEAYSTLLANMGSLSQAPGHKVEPP) of mouse Ataxin-1. Rat: 100% identity (34/34 amino acids identical). Human: 88% identity (30/34 amino acids identical).
Target	Ataxin 1
Preservatives	PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-455 was sufficient for detection of Ataxin-1 in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG2b
Clonality	Monoclonal
MW	85kDa
Uniprot ID	P54254
NCBI	NP_001186233.1
Entrez	20238
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.