

Product Datasheet

LGI1 antibody [Out of stock] (orb150840)

Description

Mouse monoclonal to LGI1 (ATTO-488).

(LGI1) was first id...

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-488

Tested

ICC, IHC, WB

Applications

Immunogen

Fusion protein amino acids 37-113 (LRRNT domain and first LRR repeat) of mouse LGI1. Rat: 100% identity (77/77 amino acids identical). Human: 98% identity (76/77 amino acids identical). ~50% identity with LGI2, LGI3 and LGI4.

Target

LGI1

Preservatives

PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-461 was sufficient for detection of LGI1 in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2a

Clonality

Monoclonal

MW

60kDa

Uniprot ID

[Q9JIA1](#)

NCBI

[NP_064674.1](#)

Entrez

[56839](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.