

Product Datasheet

QKI (pan) antibody [Out of stock] (orb150947)

Description

Mouse monoclonal to QKI (ATTO-390).

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

ATTO-390

Tested

ICC, IHC, WB

Applications
Immunogen

Fusion protein amino acids 1-341 (full-length) of human QKI-5. Mouse: 100% identity (341/341 amino acids identical). Rat: 99% identity (339/341 amino acids identical) > 90% identity with QKI-6, QKI-7 and QKI-7b.

Target

QKI

Preservatives

PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SMC-467 was sufficient for detection of Pan-QKI in 20 µg of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.

Isotype

IgG2b

Clonality

Monoclonal

MW

36-38kDa

Uniprot ID

[Q96PU8](#)

NCBI

[NP_006766.1](#)

Entrez

[9444](#)

Dilution Range

WB (1:1000)

Expiration Date

12 months from date of receipt.