

Product Datasheet

Ataxin 1 antibody (orb151092)

Description

Mouse monoclonal to Ataxin 1 (ATTO-488).

	family and contai...
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-488
Tested Applications	ICC, WB
Immunogen	Synthetic peptide amino acids 746-761 (RKRRWSAPETRKLEKS) of mouse ataxin-1. Rat: 93% identity (15/16 amino acids identical). Human: 87% identity (14/16 amino acids identical).
Target	Ataxin 1
Preservatives	PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-475 was sufficient for detection of Ataxin 1 in 20 µg of COS cells transiently transfected with Ataxin-1 lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Isotype	IgG1
Clonality	Monoclonal
MW	85kDa
Uniprot ID	P54254
NCBI	NP_001186233.1
Entrez	20238
Dilution Range	WB (1:1000)
Expiration Date	12 months from date of receipt.