

Product Datasheet

TCP1 alpha antibody (orb151177)

Description

Rat monoclonal antibody against TCP1 alpha

Species/Host	Rat
Reactivity	Bovine, Canine, Drosophila, Guinea pig, Hamster, Human, Monkey, Mouse, Other, Porcine, Rabbit, Rat, Saccharomyces, Yeast
Conjugation	PerCP
Tested Applications	FC, ICC, IP, WB
Immunogen	Recombinant Mouse TCP1 alpha protein fragment (carboxy terminal region).
Target	TCP1 alpha
Preservatives	PBS pH 7.4, 50% glycerol, 0.1% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-479 was sufficient for detection of TCP1 alpha in 20 µg of 3T3 cell lysate by colorimetric immunoblot analysis using Goat anti-rat IgG:HRP as the secondary antibody.
Isotype	IgG2a
Clonality	Monoclonal
MW	60kDa
Uniprot ID	P11983
NCBI	NP_038714.2
Entrez	21454
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.