

Product Datasheet

HSP40 antibody [Out of stock] (orb151180)

Description	Rabbit polyclonal to Hsp40 (ATTO-390). DnaJ/Hsp40 proteins have been preserved throughout evolution ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	ATTO-390
Tested Applications	ICC, IP, WB
Immunogen	Recombinant purified full-length human HSP40 (no tags)
Target	HSP40
Preservatives	Rabbit Antiserum *Storage buffer may change when conjugated
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	0.5 µg/ml of SPC-100 was sufficient for detection of HSP40 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	40kDa
Uniprot ID	P25685
NCBI	NP_006136.1
Entrez	3337
Dilution Range	WB (1:2000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.