

Product Datasheet

HSP40 antibody (orb151190)

Description

Rabbit polyclonal to Hsp40 (Biotin). DnaJ/Hsp40

evolution an...

Species/Host

Rabbit

Reactivity

Human

Conjugation

Biotin

Tested

ICC, IP, WB

Applications
Immunogen

Recombinant purified full-length human HSP40 (no tags)

Target

HSP40

Preservatives

Rabbit Antiserum *Storage buffer may change when conjugated

Concentration

N/A

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.5 µg/ml of SPC-100 was sufficient for detection of HSP40 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

40kDa

Uniprot ID
[P25685](#)
NCBI
[NP_006136.1](#)
Entrez
[3337](#)
Dilution Range

WB (1:2000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.