

## Product Datasheet

### HSP40 antibody (orb151194)

**Description**

Rabbit polyclonal to Hsp40 (PerCP). DnaJ/Hsp40

evolution and...

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

PerCP

**Tested**

ICC, IP, WB

**Applications**
**Immunogen**

Recombinant purified full-length human HSP40 (no tags)

**Target**

HSP40

**Preservatives**

Rabbit Antiserum \*Storage buffer may change when conjugated

**Concentration**

N/A

**Storage**

Conjugated antibodies should be stored at 4°C

**Note**

For research use only

**Application notes**

0.5 µg/ml of SPC-100 was sufficient for detection of HSP40 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

**Clonality**

Polyclonal

**MW**

40kDa

**Uniprot ID**
[P25685](#)
**NCBI**
[NP\\_006136.1](#)
**Entrez**
[3337](#)
**Dilution Range**

WB (1:2000), ICC/IF (1:100)

**Expiration Date**

12 months from date of receipt.