

Product Datasheet

GRP94 antibody (orb151209)

Description	Rabbit polyclonal to GRP94 (HRP). Grp94 (glucose regulated protein 94, gp96) is a constitutively exp...
Species/Host	Rabbit
Reactivity	Bovine, Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ICC, IHC, IP, WB
Immunogen	Synthetic peptide corresponding to the sequence near the C-terminus of mouse GRP94
Target	GRP94
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SPC-101 was sufficient for detection of Grp94 in 20 µg of Hela lysate by colorimetric immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	94kDa
Uniprot ID	P08113
NCBI	NP_035761.1
Entrez	22027
Dilution Range	WB (1:1000), ICC/IF (1:120), IP (1:80)
Expiration Date	12 months from date of receipt.