

Product Datasheet

HSP70 antibody (orb151232)

Description

Rabbit polyclonal to Hsp70 (ATTO-488). Hsp70

	hsps Hsp70 ...
Species/Host	Rabbit
Reactivity	Bovine, Canine, Coral, Fish, Guinea pig, Hamster, Human, Mammal, Monkey, Mouse, Plant, Porcine, Rat, Sheep
Conjugation	ATTO-488
Tested Applications	ELISA, FC, ICC, IHC, IP, WB
Immunogen	Recombinant Full length HSP70 Protein
Target	HSP70
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-103 was sufficient for detection of HSP70 in 20 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	70kDa
Uniprot ID	P0DMV8 , P0DMV9
NCBI	NP_005336.3
Entrez	3303
Dilution Range	WB (1:1000), IHC (1:100), ICC/IF (1:100), IP (1:100)
Expiration Date	12 months from date of receipt.