

Product Datasheet

HSP70 antibody (orb151240)

Description

Rabbit polyclonal to Hsp70 (APC). Hsp70 genes

	Hsp70 genes...
Species/Host	Rabbit
Reactivity	Bovine, Canine, Coral, Fish, Guinea pig, Hamster, Human, Mammal, Monkey, Mouse, Plant, Porcine, Rat, Sheep
Conjugation	APC
Tested Applications	ELISA, FC, ICC, IHC, IP, WB
Immunogen	Recombinant Full length HSP70 Protein
Target	HSP70
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-103 was sufficient for detection of HSP70 in 20 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	70kDa
Uniprot ID	PODMV8, PODMV9
NCBI	NP_005336.3
Entrez	3303
Dilution Range	WB (1:1000), IHC (1:100), ICC/IF (1:100), IP (1:100)
Expiration Date	12 months from date of receipt.