

Product Datasheet

HSP90 antibody (orb151260)

Description	Rabbit polyclonal to Hsp90 (HRP). Hsp90 is a highly conserved and essential stress protein Hsp90 is ...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, IP, WB
Immunogen	Full length protein HSP90
Target	HSP90
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:500 dilution of SPC-104 was sufficient for detection of 0.2 mg of purified HSP90 by ECL immunoblot analysis.
Clonality	Polyclonal
MW	90kDa
Uniprot ID	P08238
NCBI	NP_031381.2
Entrez	3326
Dilution Range	WB (1:500), IHC (1:100), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.