

Product Datasheet

HSP90 antibody (orb151262)

Description Rabbit polyclonal to Hsp90 (PerCP). Hsp90 is a highly conserved and essential stress protein
Hsp90 i...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation PerCP

Tested Applications ELISA, ICC, IHC, IP, WB

Immunogen Full length protein HSP90

Target HSP90

Preservatives PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration 1 mg/ml

Storage Conjugated antibodies should be stored at 4°C

Note For research use only

Application notes A 1:500 dilution of SPC-104 was sufficient for detection of 0.2 mg of purified HSP90 by ECL immunoblot analysis.

Clonality Polyclonal

MW 90kDa

Uniprot ID [P08238](#)

NCBI [NP_031381.2](#)

Entrez [3326](#)

Dilution Range WB (1:500), IHC (1:100), ICC/IF (1:100)

Expiration Date 12 months from date of receipt.