

Product Datasheet

HSP60 antibody (orb151266)

Description

Rabbit polyclonal to Hsp60 (ATTO-488). In both

In b...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Hamster, Human, Mouse, Rabbit, Rat

Conjugation

ATTO-488

Tested

ELISA, IHC, IP, WB

Applications
Immunogen

Human HSP60 produced through recombinant DNA methods in E.coli

Target

HSP60

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-105 was sufficient for detection of HSP60 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using goat anti-mouse IgG as the secondary antibody.

Clonality

Polyclonal

MW

60kDa

Uniprot ID
P10809
NCBI
NP_002147.2
Entrez
3329
Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.