

Product Datasheet

HSP60 antibody (orb151275)

Description

Rabbit polyclonal to Hsp60 (Biotin). In both

In bot...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Hamster, Human, Mouse, Rabbit, Rat

Conjugation

Biotin

Tested

ELISA, IHC, IP, WB

Applications
Immunogen

Human HSP60 produced through recombinant DNA methods in E.coli

Target

HSP60

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-105 was sufficient for detection of HSP60 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using goat anti-mouse IgG as the secondary antibody.

Clonality

Polyclonal

MW

60kDa

Uniprot ID
P10809
NCBI
NP_002147.2
Entrez
3329
Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.