

Product Datasheet

HSP27 antibody (orb151293)

Description	Rabbit polyclonal to Hsp27 (FITC). Hsp27s belong to an abundant and ubiquitous family of small heat ...
Species/Host	Rabbit
Reactivity	Canine, Fish, Human, Mouse
Conjugation	FITC
Tested Applications	FACS, FC, ICC, IHC, WB
Immunogen	Human HSP27, His tagged
Target	HSP27
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1.38 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:2000 dilution of SPC-106 was sufficient for detection of HSP27 in 20 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	27kDa
Uniprot ID	P04792
NCBI	NP_001531.1
Entrez	3315
Dilution Range	WB (1:2000), ICC/IF (1:250)
Expiration Date	12 months from date of receipt.