

Product Datasheet

GRP78 antibody (orb151313)

Description	Rabbit polyclonal to GRP78 (PerCP). GRP78 is a ubiquitously expressed, 78-kDa glucose regulated prot...
Species/Host	Rabbit
Reactivity	Bovine, Canine, Frog, Fungi, Hamster, Human, Monkey, Mouse, Plant, Rabbit, Rat
Conjugation	PerCP
Tested Applications	ICC, IHC, IP, WB
Immunogen	Rat GRP78 (Bip) synthetic peptide conjugated to KLH
Target	GRP78
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-107 was sufficient for detection of Grp78 in 10 µg of rat tissue lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	78kDa
Uniprot ID	P06761
NCBI	NP_037215.1
Entrez	25617
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.