

Product Datasheet

KDEL antibody (orb151347)

Description	Rabbit polyclonal to KDEL (PerCP). The endoplasmic reticulum is part of a protein sorting pathway, o...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	PerCP
Tested Applications	ICC, IHC, WB
Immunogen	KDEL containing peptide immunogen
Target	KDEL
Preservatives	PBS pH 7.2, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-109 was sufficient for detection of KDEL-containing proteins in 20 µg of HeLa cell lysate by ECL immunoblot analysis using goat anti-mouse IgG as the secondary.
Clonality	Polyclonal
Dilution Range	WB (1:1000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.