

Product Datasheet

HO-1 antibody [Out of stock] (orb151350)

Description

Rabbit polyclonal to HO-1 (ATTO-390).

that catalyzes Heme-oxyg...

Species/Host

Rabbit

Reactivity

Canine, Human, Mouse, Rat

Conjugation

ATTO-390

Tested

ICC, IHC, WB

Applications
Immunogen

Human heme-oxygenase (HO-1) synthetic multiple antigenic peptide

Target

HO-1

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µg/ml of SPC-112 was sufficient for detection of HO-1 in 10 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

32kDa

Uniprot ID
[P09601](#)
NCBI
[NP_002124.1](#)
Entrez
[3162](#)
Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.