

Product Datasheet

PDI antibody (orb151387)

Description

Rabbit polyclonal to PDI (ATTO-594). The three

proteins ...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Frog, Guinea pig, Hamster, Human, Invertebrate, Mouse, Mussel, Porcine, Rat, Sheep

Conjugation

ATTO-594

Tested

ICC, IHC, IP, WB

Applications

Immunogen

Rat PDI synthetic peptide (409-509 aa C-terminal) conjugated to KLH

Target

PDI

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-114 was sufficient for detection of PDI in 20 µg of HeLa cell lysate by ECL immunoblot analysis.

Clonality

Polyclonal

MW

58kDa

Uniprot ID

[P04785](#)

NCBI

[NP_001099245.2](#)

Entrez

[287164](#)

Dilution Range

WB (1:1000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.