

Product Datasheet

SOD (Mn) antibody (orb151463)

Description

Rabbit polyclonal to SOD (FITC). Superoxide

intracellula...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Rodent, Sheep

Conjugation

FITC

Tested

ELISA, IHC, IP, WB

Applications
Immunogen

Recombinant Human Mn SOD Protein

Target

SOD2

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.2 µg/ml of SPC-118 was sufficient for detection of Mn SOD in 20 µg of rat brain tissue extract by colorimetric immunoblot analysis using Goat anti-mouse IgG:AP as the secondary antibody.

Clonality

Polyclonal

MW

25kDa

Uniprot ID
P04179
NCBI
NP_000627.2
Entrez
6648
Dilution Range

WB (1:5000), IHC (1:100)

Expiration Date

12 months from date of receipt.