

Product Datasheet

Ubiquitin antibody (orb151470)

Description

Rabbit polyclonal to Ubiquitin (ATTO-488).

	eukaryoti...
Species/Host	Rabbit
Reactivity	Bovine, Canine, Drosophila, Fish, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep, Yeast
Conjugation	ATTO-488
Tested Applications	ICC, IHC, IP, WB
Immunogen	Native bovine Ubiquitin, conjugated to KLH
Target	Ubiquitin
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-119 was sufficient for detection of free ubiquitin in 15 µg of HeLa lysate by ECL immunoblot analysis using Donkey anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	10kDa
Uniprot ID	P0CG53
NCBI	NP_776558.1
Entrez	281370
Dilution Range	WB (1:1000), IHC (1:100), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.