

Product Datasheet

Ubiquitin antibody (orb151481)

Description

Rabbit polyclonal to Ubiquitin (HRP). Ubiquitin

Ubi...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Drosophila, Fish, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep, Yeast

Conjugation

HRP

Tested

ICC, IHC, IP, WB

Applications
Immunogen

Native bovine Ubiquitin, conjugated to KLH

Target

Ubiquitin

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-119 was sufficient for detection of free ubiquitin in 15 µg of HeLa lysate by ECL immunoblot analysis using Donkey anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

10kDa

Uniprot ID
P0CG53
NCBI
NP_776558.1
Entrez
281370
Dilution Range

WB (1:1000), IHC (1:100), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.