

Product Datasheet

Calreticulin antibody [Out of stock] (orb151506)

Description

Rabbit polyclonal to Calreticulin (ATTO-

highly conserved Ca...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep

Conjugation

ATTO-594

Tested

FC, ICC, IHC, IP, WB

Applications
Immunogen

Human calreticulin synthetic peptide with a cysteine residue added and the peptide conjugated to KLH

Target

Calreticulin

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-122 was sufficient for detection of Calreticulin in 20 µg of HeLa cell lysate by ECL immunoblot analysis.

Clonality

Polyclonal

MW

63kDa

Uniprot ID

P27797

NCBI

NP_004334.1

Entrez

811

Dilution Range

WB (1:1000), IHC (1:100), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.