

Product Datasheet

Calreticulin antibody (orb151517)

Description

Rabbit polyclonal to Calreticulin (PerCP).

conserved Ca2+ ...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep

Conjugation

PerCP

Tested

FC, ICC, IHC, IP, WB

Applications

Immunogen

Human calreticulin synthetic peptide with a cysteine residue added and the peptide conjugated to KLH

Target

Calreticulin

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-122 was sufficient for detection of Calreticulin in 20 µg of HeLa cell lysate by ECL immunoblot analysis.

Clonality

Polyclonal

MW

63kDa

Uniprot ID

[P27797](#)

NCBI

[NP_004334.1](#)

Entrez

[811](#)

Dilution Range

WB (1:1000), IHC (1:100), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.