

Product Datasheet

Calnexin-NT antibody (orb151564)

Description	Rabbit polyclonal antibody against Calnexin conjugated to Biotin
Species/Host	Rabbit
Reactivity	Bovine, Canine, Frog, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep
Conjugation	Biotin
Tested Applications	ICC, IHC, IP, WB
Immunogen	A 19 residue synthetic peptide based on dog calnexin and the peptide coupled to KLH
Target	Calnexin
Preservatives	Rabbit Antiserum *Storage buffer may change when conjugated
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:5000 dilution of SPC-127 was sufficient for detection of Calnexin in 20 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	90kDa
Uniprot ID	P24643
NCBI	NP_001003232.1
Entrez	403908
Dilution Range	WB (1:5000), IHC (1:100), ICC/IF (1:100), IP (1:100)
Expiration Date	12 months from date of receipt.